

Establishing CRISPR/Cas9 in *Lipomyces starkeyi*

Zoe Lau¹, Bonnie Mcneil¹, David Stuart¹


¹Department of Biochemistry, University of Alberta

Abstract

The goal of this project was to adapt the *Yarrowia lipolytica* plasmid based CRISPR/Cas9 system for usage in *Lipomyces starkeyi*. *Lipomyces starkeyi* is an oleaginous yeast, which synthesizes and stores high amounts of intracellular lipids. This specific yeast can store lipids at concentrations higher than 60% of its dry cell weight. Due to these high concentrations of lipids, *L. starkeyi* is a desired organism for the production of biofuels and other oleochemicals. However, there is a lack of knowledge and of genetic tools when trying to engineer the cells to produce these lipids for our use. The genome editing tool, CRISPR/Cas9 is efficient and simple, therefore desirable for the engineering of *L. starkeyi*. The goal was achieved by replacing the *Y. lipolytica* promoter with a *L. starkeyi* promoter, inserting guide RNA, as well as confirming cas9 protein expression.

Key words:

CRISPR, *Lipomyces starkeyi*, *Lipomyces starkeyi* CRISPR, gene editing, biofuels, *Lipomyces starkeyi* biofuels,



Establishing CRISPR/Cas9 in *Lipomyces starkeyi*

Zoe Lau, Bonnie McNeil, David Stuart

Department of Biochemistry, University of Alberta




Supported By:



Introduction

Lipomyces starkeyi is an oleaginous yeast, meaning that it synthesizes and stores high amounts of intracellular lipids. This specific yeast can store lipids at concentrations higher than 60% of its dry cell weight.¹ Due to these high concentrations of lipids, *L. starkeyi* is a desired organism for the production of biofuels and other oleochemicals. However, there is a lack of knowledge and of genetic tools when trying to engineer the cells to produce these lipids for our use. The genome editing tool, CRISPR/Cas9 is efficient and simple, therefore desirable for the engineering of *L. starkeyi*.² The goal of this project is to adapt the *Yarrowia lipolytica* plasmid based CRISPR/Cas9 system for usage in *L. starkeyi*.



Methods

Replacing the *Y. lipolytica* promoter with *L. starkeyi* *P_{PK1}*

- The PYK1p was amplified from *L. starkeyi* genomic DNA template using PCR.
- An AscI and SmaI restriction digest was done on the pCRISPR1 (*Yarrowia* /lipolytica optimized) plasmid to cut out the promoter in order to insert the new PYK1p.
- An AscI and EcoRV restriction digest was done on the PYK1p.
- The PYK1p and pCRISPR1 were ligated together, now referred to as pCRISPR1s.
- The ligated pCRISPR1s was transformed in *E. coli* and plated for colonies.
- Plasmid candidates were purified and verified by restriction digest (figure 1c).

Results

figure 1a. Strategy for constructing pCRISPR1s and insertion of guide RNA




figure 1b. pCRISPR1 AscI & SmaI restriction digest gel




figure 1c. pCRISPR1s ligated plasmid candidates





figure 1d. pCRISPR1s AvrII restriction digest gel



figure 1e. Gibson Assembly reaction pCRISPR1s




Conclusions

- pCRISPR1s vector was successfully constructed.
- PYK1p promoter was successfully inserted
- gRNA insertion requires validation
- The homology donor was successfully digested.
- Western Blot indicated expression of the cas9 protein.

Future Directions

Use CRISPR/Cas9 to knockout *Lip1* as proof of concept



Literature Cited

- (1) Lin, J.; Shen, H.; Tan, H.; Wu, S.; Hu, C.; Zhao, ZK. (2011) Lipid production by *Lipomyces starkeyi* cells in glucose solution without auxiliary nutrients. [http://dx.doi.org/10.31233/osf.io/2010.10](http://dx.doi.org/10.31233/osf.io/2010.10.10)
- (2) Shabot, S.; & Rubin, Z. (2017) Metabolic Engineering of Oleaginous Yeasts for Production of Fuels and Chemicals. <http://dx.doi.org/10.3389/fmicb.2017.02285>
- (3) Schweitz, C.; Hussain, M.; Birrell, M.; Wheadon, I. (2015) Synthetic RNA Polymerase III Promotes Facile High-Efficiency CRISPR-Cas9-Mediated Genome Editing in *Yarrowia lipolytica*. <http://dx.doi.org/10.1101/acsybio.15.00052>

Acknowledgements

- Thank you to my amazing lab team for teaching me so much this summer and giving me an enriching experience.
- Thank you to my sponsors for making this program possible for me to be a part of.
- Thank you to the WISEST team for providing this awesome program and for continuing to make positive change for women in STEM!

figure 2. Western blot indicating that the strains express the *Cas9* gene



Confirming protein expression in *L. starkeyi* strains transformed by Agrobacterium tumefaciens

- Strains were grown up to midlog phase.
- Proteins were separated via TCA, trichloroacetic acid)
- Proteins were separated by SDS-PAGE (sodium dodecyl sulfate - polyacrylamide gel electrophoresis)
- A semi dry transfer to a PVDF membrane was completed.
- The Western blot was performed with FLAG-HRP (horseradish peroxidase) antibodies and enhanced with chemiluminescence to expose on x-ray film (figure 2)

